


Results tabulation for 15679/SE PT0145: Avian Reticuloendotheliosis (REV)
Distribution Date: 23/01/2024

ELISA	ELISA kit used	Pos/Neg Cutoff value	24/7016		24/7017		24/7018		24/7019	
			SP Ratio	Interpretation (+ve/-ve/IC)	SP Ratio	Interpretation (+ve/-ve/IC)	SP Ratio	Interpretation (+ve/-ve/IC)	SP Ratio	Interpretation (+ve/-ve/IC)
Intended			1/8 dilution		Negative		1/16 dilution		1/4 dilution	
				+ve		-ve		+ve		+ve
1143	BioChek FS 9264	>=0.5	2,7222	+ve	0,219	-ve	1,946	+ve	3,601	+ve
1318	REV IDEXX : FW276	0.501	1.659	+ve	0.006	-ve	0.906	+ve	2.702	+ve
1923	IDEXX FW276	<=0.50 Negative, >0.50 Positive	2.42	+ve	0.07	-ve	1.89	+ve	5.05	+ve

ELISA	24/7020	
	SP Ratio	Interpretation (+ve/-ve/IC)
Intended	1/10 dilution	
		+ve
1143	1,891	+ve
1318	0.812	+ve
1923	1.46	+ve

Testing carried out to determine the intended results of the samples using IDEXX ELISA, was subcontracted to an Animal and Plant Health Agency (APHA) laboratory.

Comments provided by Julieanne Caldwell, APHA Lasswade

All three participating laboratories obtained the intended results for all samples.



0004

Version: 21/12/2017

Paula Johnson
Scheme Coordinator
Animal and Plant Health Agency
Quality Assurance Unit
The Elms, College Road
Sutton Bonington
Loughborough
Leicestershire
LE12 5RB
UK
Authorised: 05/03/2024